## Can Your Force a Graph?

Materials: Force Plate, Motion Detector, Go!Link, Logger Lite, Computer, Fathom

**Groups:** Teachers need to be in groups of 3-4. Roles include: experimenter to apply force to force plate, person to set up and turn on experiment, and recorder (to have activity open on their computer).

# **Mathematical goals:**

- For a function that models a relationship between two quantities, interpret key features of graphs and tables in terms of the quantities, and sketch graphs showing key features given a verbal description of the relationship. *Key features include: intercepts; intervals where the function is increasing, decreasing, positive, or negative; relative maximums and minimums; symmetries; end behavior; and periodicity.* (F-IF.4.)
- Relate the domain of a function to its graph and, where applicable, to the quantitative relationship it describes. For example, if the function h(n) gives the number of person-hours it takes to assemble n engines in a factory, then the positive integers would be an appropriate domain for the function. (F.IF.5.)

## **Technological goals:**

- Teachers will be able to match forces applied to a force plate to a graph
- Teachers will be able to import their data into Fathom
- Teachers will be able to create tables and graphs in Fathom based on their data

#### Lesson Outline:

#### **Background:**

Weight is a force based on the mass of an object times it's acceleration due to gravity. This comes directly from Newton's Second Law of Motion: *F=ma*. This is why objects weigh less on the moon. In this activity you'll try to match the graph of forces over time using a force plate.

## Set up the experiments in Logger Lite.

- 1) Open a new Logger Lite Window
- 2) Connect the force plate to the Go!Link cable and connect to your computers USB port.
- 3) Connect your Motion Detector to your computer. (Logger Lite only activates the Graph Match feature when the Motion Detector is attached, but we can trick it.)
- 4) Once you see the Match icon on the menu bar you can unplug the Motion Sensor.

4) Click on Match on the top menu bar. It will give you a graph to match in an experiment.

# **Make Predictions and Conjectures:**

Predict what you need to do to create:

- a) A line that is increasing
- b) A line that is decreasing
- c) A line that is constant (horizontal)
- d) A graph with a given y-intercept
- e) A graph with a negative y-intercept
- f) A graph that begins at zero
- g) A graph with a steep positive slope
- h) A graph with a gradual negative slope

### Matching Graphs and Making Sense of the Data

Discuss with your group the following questions. Record your answers.

- 1) Try to match various graphs to your data. What part is most challenging?
- 2) Were your predictions correct?
- 3) What did the different parts of your graph represent?
- 4) Are there any graphs you could draw that would be impossible to model with the force plate?
- 5) How could you use this activity in your Algebra 1 classes?

#### Conclusion

In this activity you tried to match graphs by applying forces to a force plate. You experienced what it means to increase a force, decrease a force, and how to keep it constant. You also learned to model steeper and gradual slopes of lines.

#### For Further Exploration:

Draw graphs on paper and try to see if your group member can model them. Why might some of these graphs be impossible?